



Informed
conversations
at the
point of care

HbA1c results in 5 minutes

- Improve patients' in-office experience
- Increase efficiency
- Help drive better health outcomes for diabetes patients

A1CNow⁺



A1CNow⁺ Test System Key Features

- Lab-quality accuracy
- Portable
- Easy to use
- Fingerstick (5 μ L) blood sample
- NGSP-certified
- IFCC-traceable
- CLIA-waived

A1CNow⁺ Controls™ and Connectivity



Ensure Accurate A1C Results

A1CNow⁺ Controls™
REF 737

- Ready-to-use control materials
- Two-level control kit
- Dropper vials



Room Temperature Storage

7 Days
once thawed and opened
18° to 25°C (64° to 77° F)

Refrigerated Storage

8 Months
once thawed and opened
2° to 8°C (36° to 46° F)


Long-Term Frozen Storage

3 Years
or until expiration date
-30° to -10°C (-22° to 14° F)

PTS Connect Blu-Dock™

(Sold separately)

REF 4715



- Wirelessly transfers results to third-party data management systems¹
- Bluetooth enabled
- Seamlessly connects with a PTS Connect Printer

PTS Connect Printer

(Sold separately)

REF 2750



- Outputs results on paper or self-adhesive labels
- Can connect to A1CNow⁺ through Blu-Dock

Available in a 10-count or 20-count test kit



A1CNow⁺ system 10-count test kit
REF 3024



A1CNow⁺ system 20-count test kit
REF 3021

¹ May require additional development.

For more information or to place an order:



AmericanHospitalSupply.com
Email: salesinquiries@americanhospitalSupply.com